Porcelain Liner M

EN

INSTRUCTIONS

IMPORTANT :
READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.
KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Adhesive Primer for Porcelain Bonding

Porcelain Liner M

FOR DENTIST OR DENTAL TECHNICIAN USE ONLY

1.What is Porcelain Liner M?

Porcelain Liner M is a silane coupling agent for porcelain surface. Application of the product on a porcelain surface before bonding with a dental adhesive resin cement such as Super-Bond C&B will improve the bond strength and durability.

2.Contents of the kit

- ①Liquid A ----- 5mL
- ②Liquid B ----- 5mL
- ③Plastic Dispensing Dish ----- 15
- ④Sponge (L-S) ----- 1
- Attached document
- Instructions

3.Precautions : Read all instructions thoroughly before use.

3.1 Safety Please keep the following precautions for safe use.

- ①**Applications**
Use Porcelain Liner M only for the applications recommended in this leaflet.
- ②**Past history of sensitivity**
Porcelain Liner M should not be used by clinicians / technicians or on patients who are methacrylic monomer-sensitive.
- ③**Symptomatic irritation**
Cease using Porcelain Liner M immediately, if signs of irritation such as rashes appear,and see a doctor.
- ④**Avoid contact**
Avoid contact with soft tissue,skin or eyes. A rubber dam is recommended for intraoral use. Contaminated skin or mucosa should be wiped off immediately with alcohol and then thoroughly rinsed with running water. If Porcelain Liner M enters the eye, immediately rinse thoroughly with running water. The patient should be examined by an ophthalmologist. Be sure to wear dental gloves.
- ⑤**Give care to flammability**
Liquid A and Liquid B are flammable. Do not store where they may be exposed to open flame.

3.2 Storage

- ①**Storage conditions**
Store the Liquid A and Liquid B in a cool,dark location. High temperature,high humidity and direct sunlight will shorten their shelf-life.
- ②**Volatility**
Liquid A and Liquid B are highly volatile. Recap the bottles immediately.
- ③**Contamination**
Do not mix the bottle caps. Plastic Dispensing Dish is disposable. Do not reuse, and avoid using with other materials.

As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.

4. How to use Porcelain Liner M

- ①**Pretreatment of the porcelain surface**
Remove contaminants and stains fromthe porcelain surface.

Depending on the case at hand,clean by use of abrasive wheels and points or a prophy cup with fluoride-free,oil-free pumice, then rinse thoroughly with water and dry.

*Applying the Red Activator will help remove surface contaminants.

*Ultrasonic cleaner is also effective.

②Preparation of the mixture of Liquid A and Liquid B

Dispense Liquid A and Liquid B to a well of a Plastic Dispensing Dish and mix with an application sponge.

[Mixing Ratio]

Liquid A / Liquid B = 1 / 1

*The mixture becomes usable immediately.

*The mixture is volatile. Use quickly after mixing.

③Application of Porcelain Liner M

Using an application sponge in the kit, apply the mixture to the porcelain surface to be bonded.

*Blow lightly if the liquid remains in a dent.

*Application of the mixture to tooth or metal surface will not affect the adhesion.

*If possible, after applying Porcelain Liner M, heat the surface with a hair dryer or a warm-air dryer to enhance the bond strength.

[Treating Condition]

At 80-120°C for 2-3 min.

④Bonding Operation

Apply the adhesive resin cement following the instructions.

*Carefully avoid contamination of the applied surface, and apply the resin cement quickly.

5.Clinical Applications

- ①Repair of porcelain
 - Repair of a fractured porcelain tooth
 - Repair with composite resin of a fractured porcelain veneered crown
 - Adhesion of all-ceramic crown or inlay
 - Adhesion of porcelain veneer

- ②Adhesion of orthodontic bracket to porcelain tooth

- ③Adhesion of ceramic bracket to tooth surface

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Made in Japan by  **SUN MEDICAL CO.,LTD.** 571-2, Furutaka-cho, Moriyama, Shiga, 524-0044, Japan
EC-Representative : **MITSUI CHEMICALS EUROPE GmbH**